

RELM beta antibody

Cat. No. GTX30435

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent

Not tested in other applications.

Calculated MW 11 kDa. ([Note](#))

Product Note The sample used was secreted mouse RELM Beta from the media of transfected 293T cells. It should be pointed out that mouse RELM beta is a secreted protein. This might preclude finding a signal in tissue extract, even though secreted protein.

Properties

Form	Liquid
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Reacts with the residues 32-46 of the mouse RELM beta protein
Purification	IgG
Conjugation	Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx30435 WB Image

Western blot analysis of secreted mouse RELM beta from the media of transfected 293T cells with anti-mouse RELM beta at a 1:500 dilution



For full product information, images and publications, please visit our [website](https://www.genetex.com).