

ABCG8 antibody

Cat. No. GTX30466

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

References (1)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500

Not tested in other applications.

Calculated MW 76 kDa. ([Note](#))

Product Note We do not recommend use of this product for Mouse samples.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic peptide from the N-terminal region of human ABCG8 protein. [UniProt# Q9H221]
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX30466 WB Image

WB analysis of mouse liver microsomes from abcg5/abcg8 knockout mice, wildtype human abcg5/abcg8, and transgenic human abcg5/abcg8 using GTX30466 ABCG8 antibody.



GTX30466 WB Image

The data was published in the journal PLoS One in 2016. [PMID: 27023132](https://pubmed.ncbi.nlm.nih.gov/27023132/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).