


# Antigen retrieval solution

**Cat No. GTX30709****Application** IHC-P**Package**  
100 ml**PRODUCT****Summary** Antigen retrieval - heat mediated pH10**APPLICATION****Application Note**

IHC-P: 1X buffer solution is intended for heat-induced antigen retrieval in IHC. See protocol

**PROPERTIES****Form** Liquid**Buffer** 10X buffer solution**Storage** Keep as concentrated solution, aliquot and store at 4°C. DO NOT FREEZE.**Note** For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.**DATA IMAGES****GTX30709 IHC-P Image**For full product information, images and publications, please visit our [website](#).