

## Antigen retrieval solution

Cat. No. GTX30709

Application	IHC-P
-------------	-------

Package	100 ml
---------	--------

## PRODUCT

Summary	Antigen retrieval - heat mediated pH10
---------	--

## APPLICATION

## Application Note

IHC-P: 1X buffer solution is intended for heat-induced antigen retrieval in IHC. See protocol

## PROPERTIES

Form	Liquid
------	--------

Buffer	10x buffer solution
--------	---------------------

Preservative	No preservative
--------------	-----------------


Storage	Keep as concentrated solution, aliquot and store at 4°C. DO NOT FREEZE.
---------	---

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



GTX30709 IHC-P Image



For full product information, images and publications, please visit our [website](#).