

Antigen retrieval solution

Cat. No. GTX30709

Application IHC-P

Package
100 ml

PRODUCT

Summary Antigen retrieval - heat mediated pH10

APPLICATION

Application Note

IHC-P: 1X buffer solution is intended for heat-induced antigen retriever in IHC. See protocol

PROPERTIES

Form Liquid

Buffer 10x buffer solution

Preservative No preservative


Storage Keep as concentrated solution, aliquot and store at 4°C. DO NOT FREEZE.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX30709 IHC-P Image



For full product information, images and publications, please visit our [website](#).