

Target Unmasking Fluid

Cat. No. GTX30718

Application IHC-P

Package
50 ml

PRODUCT

Summary Target Unmasking Fluid

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	Assay dependent

Not tested in other applications.

PROPERTIES

Form Liquid

Buffer Antigen unmasking fluid is supplied as 3x concentrated fluid. Crystals in the fluid will dissolve during heating before use.

Preservative No preservative

Storage Store at room temperature.

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).