

## TNP-BSA

Cat. No. GTX30726-pro

Package

1 ml

## PRODUCT

**Summary** 2, 4, 6 Trinitrophenyl (TNP) is conjugated to Bovine Serum Albumin (BSA) through lysine residues by amide bonds.

## Properties

**Form** Liquid

**Buffer** 10mM PO<sub>4</sub>, 200mM NaCl, 0.5mM EDTA

**Preservative** 0.05% Sodium azide

**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. Protect from light.

**Concentration** 5 mg/ml (Please refer to the vial label for the specific concentration.)

**Conjugation** BSA

**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).