

GPR120 peptide

Cat. No. GTX30752

Applications Neutralizing /Inhibition

Package
25 µg

Applications

Application Note

The peptide is used for blocking the antibody activity of GPCR GPR120 antibody (GeneTex Cat.# GTX71279).

Product Note

It usually blocks the antibody activity completely in immunohistochemistry analysis by pre-incubating the antibody with a 1-10X molar volume of peptide for 60 minutes at room temperature. Optimization of the molar ratio of antigen to antibody may be required.

Properties

Form Liquid

Buffer PBS

Preservative No preservatives

Storage Keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Purification =95% by HPLC

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).