

TYK2 antibody

Cat. No. GTX30812

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:2,000
ICC/IF	Assay dependent

Not tested in other applications.

Calculated MW 134 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	synthetic peptide corresponding to amino acids 1046-1060 of human TYK2.
Purification	IgG fraction
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX30812 ICC/IF Image

ICC/IF analysis of HUVEC cells using GTX30812 TYK2 antibody at 1:100 (red) with Hoescht 33342 (blue).



GTX30812 WB Image

WB analysis of HEK-293T cells overexpressing human TYK2 using GTX30812 TYK2 antibody at 1:1000.



For full product information, images and publications, please visit our [website](#).