

SLUG antibody

Cat. No. GTX30813

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

References (2)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5-2µg/ml
IHC-P	Assay dependent

Note : This antibody detects a band of approximately 30kDa in 293 cell lysate.

Not tested in other applications.

Calculated MW 30 kDa. ([Note](#))

Product Note This antibody detects Slug and is not cross-reactive with Snail.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	13 amino acid peptide sequence near the centre of human Slug.
Purification	Purified by affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES

GTx30813 IHC-P Image

IHC-P analysis of human kidney tissue using GTx30813 SLUG antibody.

Dilution : 20 µg/ml


GTx30813 WB Image

WB analysis of HEK293 cell lysate using GTx30813 SLUG antibody at 0.5(A), and 1(B) and 2(C) µg/ml.


GTx30813 WB Image

WB analysis of (A) human and (B) mouse kidney tissue lysate using GTx30813 SLUG antibody.

Dilution : 0.5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).