

## LAMP2 blocking peptide

Cat. No. GTX30881-PEP

Application	Neutralizing/Inhibition
-------------	-------------------------

Package	50 µg
---------	-------

## PRODUCT

Summary	LAMP2 peptide is used for blocking the activity of LAMP2 antibody (catalog number GTX30881).
---------	--

## APPLICATION

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

Neutralizing/Inhibition	Assay dependent
-------------------------	-----------------

Not tested in other applications.

## PROPERTIES

Form	Liquid
------	--------

Buffer	PBS, 0.1% BSA
--------	---------------

Preservative	0.02% Sodium azide
--------------	--------------------

Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
---------	--

Concentration	200 µg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Conjugation	Unconjugated
-------------	--------------

Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>
------	---



For full product information, images and publications, please visit our [website](#).