

Ly-6A/E antibody [D7]

Cat. No. GTX30898

Host	Rat
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, IHC-Fr, FCM, IP
Reactivity	Mouse

References (1)

Package

250 µg

Applications

Application Note

For FACS: Use at an assay dependent dilution. For IHC-Fr: Use at an assay dependent dilution. For IP: Use at an assay dependent dilution. For WB: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW 14 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	Borate buffered saline
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	IL-2-dependent mouse T-cell line CTL-L
Purification	Purified IgG
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).