

PLCB1 antibody [O.T.171]

Cat. No. GTX30899

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB, IP
Reactivity	Human, Mouse, Rat, Bovine

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	5µg/ml
IP	5µg

Not tested in other applications.

Calculated MW 139 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	100mM Tris-Glycine, 150mM NaCl
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified soluble PLC (phospholipase C) b-1 from bovine brain.
Purification	Protein G purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX30899 IP Image

IP analysis of bovine brain membrane preparation using GTX30899 PLCB1 antibody [0.T.171]. Arrows indicate PLC b-1 (150kD), heavy (H) and light (L) IgG chains.

IP reaction : 5µg antibody/mg lysate

Dilution : 5µg/ml



For full product information, images and publications, please visit our [website](#).