

Urease protein (Biotin)

Cat. No. GTX30905-pro

Applications Functional Assay

Package
1 ml

Applications

Application Note

Before conjugation, the activity was ≥ 45 Units/mg.

Product Note Reacts with immobilized avidin > 95%

Properties

Form Liquid

Buffer 50mM PO₄, 150mM NaCl, 40% Glycerol, 1mM EDTA

Preservative 0.05% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 2.5 mg/ml (Please refer to the vial label for the specific concentration.)

Expression System Jack beans (*Canavalia ensiformis*)

Conjugation Biotin

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).