

Urease protein (Biotin)

Cat. No. GTX30905-pro

Applications

Functional Assay

Package

1 ml

Applications

Application Note

Before conjugation, the activity was ≥ 45 Units/mg.

Product Note

Reacts with immobilized avidin > 95%

Properties

Form	Liquid
Buffer	50mM PO ₄ , 150mM NaCl, 40% Glycerol, 1mM EDTA
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	2.5 mg/ml (Please refer to the vial label for the specific concentration.)
Expression System	Jack beans (Canavalia ensiformis)
Conjugation	Biotin
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).