

## Fluoroshield™ with DAPI

**Cat. No. GTX30920****Applications** ICC/IF, IHC-Fr, IHC

References ( 58 )

Package

20 ml

## PRODUCT

**Summary** Fluoroshield™ with DAPI

## Applications

**Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

**Suggested dilution** Recommended dilution

ICC/IF Assay dependent

IHC-Fr Assay dependent

IHC Assay dependent

Not tested in other applications.


## Properties

**Form** Liquid**Preservative** 0.09% Sodium azide**Storage** Store at 2-8°C, Protect from light, DO NOT FREEZE.**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES

**GTX30920 ICC/IF Image**

Fluoroshield™ with DAPI is an aqueous mounting medium for preserving fluorescence of tissue and cell smears. This mounting medium is fortified with DAPI which is a counterstain for DNA.

For full product information, images and publications, please visit our [website](#).