

## Fluoroshield™ with DAPI

Cat. No. GTX30920

Applications ICC/IF, IHC-Fr, IHC

References ( 50 )

Package

20 ml

## PRODUCT

Summary Fluoroshield™ with DAPI

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

ICC/IF	Assay dependent
--------	-----------------

IHC-Fr	Assay dependent
--------	-----------------

IHC	Assay dependent
-----	-----------------

Not tested in other applications.

## Properties

Form	Liquid
------	--------


Storage	Store at 2-8°C, Protect from light, DO NOT FREEZE.
---------	--

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES



## GTX30920 ICC/IF Image

Fluoroshield™ with DAPI is an aqueous mounting medium for preserving fluorescence of tissue and cell smears. This mounting medium is fortified with DAPI which is a counterstain for DNA.



For full product information, images and publications, please visit our [website](#).