

Fluoroshield™ with DAPI

Cat. No. GTX30920

Applications ICC/IF, IHC-Fr, IHC References (51)
Package
20 ml

PRODUCT

Summary Fluoroshield™ with DAPI

Applications

Application Note


*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
IHC-Fr	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Properties	
Form	Liquid
Storage	Store at 2-8°C, Protect from light, DO NOT FREEZE.
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX30920 ICC/IF Image

Fluoroshield $^{\text{TM}}$ with DAPI is an aqueous mounting medium for preserving fluorescence of tissue and cell smears. This mounting medium is fortified with DAPI which is a counterstain for DNA.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 12 / 02 Page 1 of 1

€ 886-3-6208988 📻 886-3-6208989 🐷 infoasia@genetex.com