

# ImmunoFluoroMount™

**Cat. No. GTX30928**

**Applications** IHC

References ( 6 )  
 Package  
 30 ml

## PRODUCT

### Summary

Refractive Index:  $1.400 \pm 0.002$   
 (These numbers apply to these mounting mediums in solution. Refractive indexes change when the water solvent evaporates and mounting media dries on slides. We do not have the means to measure the refractive indexes of dry mounting mediums; however we expect the numbers to go higher when dried. The refractive index of water is 1.3330)

## Applications

### Application Note

Apply 2-3 drops of mounting medium and apply coverslip carefully without getting air bubbles.

## Properties

**Form** Liquid

**Storage** Store at 2-8°C, Protect from light, DO NOT FREEZE.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.  
 Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).