

Universal Antibody Dilution Buffer™ (ready to use)

Cat. No. GTX30932

Application WB, ICC/IF, ELISA, IHC

Reference (1)

Package

100 ml

PRODUCT

Summary Universal Antibody Dilution Buffer™ (ready to use)

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
ELISA	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Product Note

This buffer does not contain any mammalian proteins, phosphate, sodium azide, or mercury preservative and can be used for dilution of all antibodies, including peroxidase and antibodies to phosphoproteins. Not suitable for dilution of antibodies to S100 proteins.

PROPERTIES

Form Liquid

Storage Store at 2-8°C.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).