

# Universal Antibody Dilution Buffer™ (ready to use with BSA, Immunoglobulin free)

**Cat. No. GTX30933**

**Applications** WB, ICC/IF, ELISA, IHC

**Package**

100 ml

## PRODUCT

**Summary** Universal Antibody Dilution Buffer™, ready to use with BSA, Immunoglobulin free

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
ELISA	Assay dependent
IHC	Assay dependent

Not tested in other applications.

## Properties

**Form** Liquid

**Preservative** No preservatives

**Storage** Store at 2-8°C, DO NOT FREEZE. The diluted antibodies can be stored at 2-8°C for several weeks in this buffer.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).