

TMBM (For Membranes)

Cat. No. GTX30946

Applications	WB
--------------	----

Package	50 ml
---------	-------

PRODUCT

Summary	TMBM (For Membranes)
---------	----------------------

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

WB	Assay dependent
----	-----------------

Not tested in other applications.

Properties

Storage	Store at 2-8°C.
---------	-----------------

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).