

TMBM (For Membranes)

Cat. No. GTX30946

Application	WB
--------------------	----

Package	50 ml
----------------	-------

PRODUCT

Summary	TMBM (For Membranes)
----------------	----------------------

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

WB	Assay dependent
----	-----------------

Not tested in other applications.

PROPERTIES

Storage	Store at 2-8°C.
----------------	-----------------

Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>
-------------	---



For full product information, images and publications, please visit our [website](#).