

TMB Substrate

Cat. No. GTX30947

Application IHC

Package
50 ml

PRODUCT

Summary TMB-IHC (IHC,ICC,In Situ)

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC	Assay dependent

Not tested in other applications.

PROPERTIES

Storage Store at 2-8°C.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).