

Tween 20 (50%)

Cat. No. GTX30962

Application WB, ICC/IF, ELISA, IHC

Reference (1)

Package

50 ml

PRODUCT

Summary Tween 20 (50%)

APPLICATION

Application Note

For diluting antibodies use 1:1000 in buffers (1 ml into 999 ml).

Washing buffer use 1:500 (1 ml into 499 ml).

The optimum dilution should be determined by the individual lab.

PROPERTIES

Buffer This detergent is supplied as 50% in water

Preservative 0.05 % preservative, not contain azide or mercury preservative

Storage Store at room temperature.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).