

BSA (IgG and protease free)

Cat. No. GTX30964

Application WB, ICC/IF, ELISA, IHC

Reference (1)

Package

10 g

PRODUCT

Summary BSA (IgG and protease free)

APPLICATION

Application Note

Reconstitute in desired buffer to the concentration or % of choice.

PROPERTIES

Form Lyophilized powder

Storage Store at 2-8°C.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).