

# Endogenous peroxidase blocking kit

**Cat. No. GTX30967****Applications** ICC/IF, IHC

References ( 4 )

Package

15 ml

**PRODUCT****Summary** Endogenous peroxidase blocking kit**Applications****Application Note**

We recommend that endogenous peroxidase activity should be blocked following deparaffinization but prior to antigen retrieval.

**Properties****Storage** Store at 2-8°C, DO NOT FREEZE.

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

**Note**

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

**DATA IMAGES****GTX30967 Image**

For full product information, images and publications, please visit our [website](#).