

Endogenous peroxidase blocking kit

Cat. No. GTX30967
Applications ICC/IF, IHC

References (4)

Package

15 ml

PRODUCT

Summary Endogenous peroxidase blocking kit

Applications

Application Note

We recommend that endogenous peroxidase activity should be blocked following deparaffinization but prior to antigen retrieval.

Properties

Storage Store at 2-8°C, DO NOT FREEZE.

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES


GTX30967 Image


For full product information, images and publications, please visit our [website](#).