

# Endogenous peroxidase blocking kit

**Cat. No. GTX30967**
**Applications** ICC/IF, IHC

**References** ( 4 )

**Package**

15 ml

## PRODUCT

**Summary** Endogenous peroxidase blocking kit

## Applications

## Application Note

We recommend that endogenous peroxidase activity should be blocked following deparaffinization but prior to antigen retrieval.

## Properties

**Storage** Store at 2-8°C, DO NOT FREEZE.

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

## DATA IMAGES


**GTX30967 Image**


For full product information, images and publications, please visit our [website](#).