

# Blocking Solution with Bovine Serum (ready to use)

**Cat. No. GTX30970**

<b>Application</b>	WB, ICC/IF, ELISA, Blocking, IHC
<b>Species</b>	Bovine

**Package**  
50 ml

## PRODUCT

**Summary** Immuno Protein Blocking Solution with Bovine Serum (ready to use)

## APPLICATION

### Application Note

This Immuno Protein Blocking solutions can be used in with all Immunoassays IHC, ICC, WB, ELISA) to reduce non-specific background and conserved sequence. It is recommended that this Immuno Blocking solution should be used from the same host animal where the labeled secondary antibody is produced; e.g. if the labeled secondary antibody is made in Chicken, the blocking should be one with Chicken serum.

## PROPERTIES

**Form** Liquid

**Storage** Store at 2-8°C.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).