

Blocking Solution with Donkey Serum (ready to use)

Cat. No. GTX30972

Applications	WB, ICC/IF, ELISA, Blocking, IHC
Species	Donkey

References (11)
Package
50 ml

PRODUCT

Summary Immuno Protein Blocking Solution with Donkey Serum (ready to use)

Applications

Application Note

The serum used to make this blocking solution is obtained from pools of healthy adult animals. The serum is heat inactivated at 56°C for 60 minutes and dilapidated and filtered. The serum is tested with appropriate biotinylated antibody for cross-reactivity. Immuno protein blocking solutions is ready to use and do not contain azide or phosphate buffers.

Properties

Form Liquid

Storage Store at 2-8°C.

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).