

# Blocking Solution with Goat Serum (ready to use)

**Cat. No. GTX30973**

<b>Application</b>	WB, ICC/IF, ELISA, Blocking, IHC
<b>Species</b>	Goat

Reference ( 3 )  
Package  
50 ml

## PRODUCT

**Summary** Immuno Protein Blocking Solution with Goat Serum (ready to use)

## APPLICATION

### Application Note

The serum used to make this blocking solution is obtained from pools of healthy adult animals. The serum is heat inactivated at 56°C for 60 minutes and diluted and filtered. The serum is tested with appropriate biotinylated antibody for cross-reactivity. Immuno protein blocking solutions is ready to use and do not contain azide or phosphate buffers.

## PROPERTIES

**Form** Liquid

**Storage** Store at 2-8°C.

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).