

RANTES antibody

Cat. No. GTX31209

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA, IHC
Reactivity	Rat

Package
100 µg

Applications

Application Note

We recommend the following starting dilutions: For WB: Use at a concentration of 0.1-0.2 µg/ml. For ELISA: Use at a concentration of 0.5-2.0 µg/ml. For IHC: Use an assay dependent dilution. Optimal dilutions should be determined experimentally by the end user.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	E.coli derived Recombinant Rat RANTES (CCL5)
Purification	Affinity Purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX31209 IHC Image

IHC analysis of colchicine injected rat brain (somatosensory cortex) tissue using RANTES antibody at a concentration of 1.0 mg/ml. This was followed by a peroxidase conjugated secondary antibody and then a fluorescein Tyramide Signal Amplification (TSA) reagent.



For full product information, images and publications, please visit our [website](#).