

DFFB antibody

Cat. No. GTX31257

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	2 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 39 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	DFFB antibody was raised against a peptide corresponding to 17 amino acids near the center of murine DFFB. The immunogen is located within amino acids 190 - 240 of DFFB.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

DATA IMAGES

GTX31257 IHC-P Image


IHC-P analysis of mouse kidney tissue using GTX31257 DFFB antibody.
Working concentration : 2 µg/ml


GTX31257 IHC-P Image

IHC-P analysis of mouse kidney tissue using GTX31257 DFFB antibody.
Working concentration : 5 µg/ml
Green : Primary antibody
Blue : DAPI
Red : Actin


GTX31257 IHC-P Image

IHC-P analysis of mouse kidney tissue using GTX31257 DFFB antibody.
Working concentration : 10 µg/ml


GTX31257 WB Image

WB analysis of mouse kidney tissue lysate using GTX31257 DFFB antibody.
Working concentration : 2 µg/ml



For full product information, images and publications, please visit our [website](#).