

IRAK3 antibody

Cat. No. GTX31265

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
ICC/IF	Assay dependent
IHC-P	2 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 68 kDa. ([Note](#))

Product Note IRAK-M antibody is predicted to have no cross reactivity to other members of the IRAK family.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	IRAK3 antibody was raised against a peptide corresponding to 16 amino acids near the carboxy terminus of human IRAK3. The immunogen is located within the last 50 amino acids of IRAK3.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX31265 IHC-P Image

IHC-P analysis of human spleen tissue using GTX31265 IRAK3 antibody.


Working concentration : 5 µg/ml



GTX31265 WB Image

WB analysis of (M) mouse spleen and (R) rat liver tissue lysates using GTX31265 IRAK3 antibody.

Working concentration : 1 µg/ml



GTX31265 IHC-P Image

IHC-P analysis of rat liver tissue using GTX31265 IRAK3 antibody.

Working concentration : 2 µg/ml



For full product information, images and publications, please visit our [website](#).