

CARD14 antibody

Cat. No. GTX31281

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Product Note

CARMA2 antibody is human specific. At least three isoforms of CARMA2 are known to exist; this antibody will only detect isoform 1. CARMA2 antibody is predicted not to cross-react with other CARMA proteins.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	CARD14 antibody was raised against a 14 amino acid synthetic peptide near the carboxy terminus of human CARD14. The immunogen is located within the last 50 amino acids of CARD14.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX31281 IHC-P Image

IHC-P analysis of human placenta tissue using GTX31281 CARD14 antibody.
Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](#).