

## TSC2 blocking peptide

Cat. No. GTX31287-PEP

## Applications

Neutralizing/Inhibition

## Package

50 µg

## PRODUCT

## Summary

TSC2 peptide is used for blocking the activity of TSC2 antibody (catalog number GTX31287).

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

## Suggested dilution

## Recommended dilution

Neutralizing/Inhibition

Assay dependent

Not tested in other applications.

## Properties

## Form

Liquid

## Buffer

PBS, 0.1% BSA

## Preservative

0.02% Sodium azide

## Storage

Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

## Concentration

200 µg/ml (Please refer to the vial label for the specific concentration.)

## Conjugation

Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).