

## CRMP1 antibody

Cat. No. GTX31288

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

**Package**  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 - 2 µg/mL
IHC-P	2 µg/mL
ELISA	Assay dependent

Not tested in other applications.

**Calculated MW** 62 kDa. ( [Note](#) )

## Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	CRMP1 antibody was raised against a 15 amino acid synthetic peptide from near the carboxy terminus of human CRMP1. The immunogen is located within the last 50 amino acids of CRMP1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).


Date 2026 / 01 / 10 Page 1 of 2

**DATA IMAGES**

**GTX31288 IHC-P Image**

IHC-P analysis of human brain tissue using GTX31288 CRMP1 antibody.  
Working concentration : 2 µg/ml


**GTX31288 IHC-P Image**

IHC-P analysis of human brain tissue using GTX31288 CRMP1 antibody.  
Working concentration : 10 µg/ml


**GTX31288 WB Image**

WB analysis of human brain tissue lysate using GTX31288 CRMP1 antibody.  
Working concentration : (A) 1 and (B) 2 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).