

BUB3 antibody

Cat. No. GTX31314

| | |
|---------------------|-------------------|
| Host | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Applications | WB, ICC/IF, ELISA |
| Reactivity | Human, Mouse |

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB | 0.5 - 1 µg/mL |
| ICC/IF | 10 µg/mL |
| ELISA | Assay dependent |

Not tested in other applications.

Calculated MW 37 kDa. ([Note](#))

Properties

| | |
|----------------------|--|
| Form | Liquid |
| Buffer | PBS |
| Preservative | 0.02% Sodium azide |
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.) |
| Immunogen | Bub3 antibody was raised against a 16 amino acid synthetic peptide from near the carboxy terminus of human bub3. The immunogen is located within the last 50 amino acids of Bub3. |
| Purification | Purified by antigen-affinity chromatography |
| Conjugation | Unconjugated |

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES



GTX31314 WB Image

WB analysis of Jurkat cell lysate using GTX31314 BUB3 antibody.


Working concentration : (A) 0.5 and (B) 1 µg/ml



GTX31314 ICC/IF Image

ICC/IF analysis of Jurkat cells using GTX31314 BUB3 antibody.

Working concentration : 20 µg/ml



GTX31314 ICC/IF Image

ICC/IF analysis of Jurkat cells using GTX31314 BUB3 antibody.

Working concentration : 10 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).