

## ATG10 antibody

Cat. No. GTX31329

|              |                   |
|--------------|-------------------|
| Host         | Rabbit            |
| Clonality    | Polyclonal        |
| Isotype      | IgG               |
| Applications | WB, ELISA         |
| Reactivity   | Human, Mouse, Rat |

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 0.5 - 2 µg/mL        |
| ELISA              | Assay dependent      |

Not tested in other applications.

Calculated MW 25 kDa. ( [Note](#) )

## Properties

|               |                                                                                                                                                                                      |
|---------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                               |
| Buffer        | PBS                                                                                                                                                                                  |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                   |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.                                                                                              |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                             |
| Immunogen     | ATG10 antibody was raised against a 15 amino acid synthetic peptide from near the carboxy terminus of human ATG10. The immunogen is located within the last 50 amino acids of ATG10. |
| Purification  | Purified by antigen-affinity chromatography                                                                                                                                          |
| Conjugation   | Unconjugated                                                                                                                                                                         |

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX31329 WB Image

WB analysis of SK-N-SH cell lysate using GTX31329 ATG10 antibody.

Working concentration : (A) 0.5, (B) 1, and (C) 2  $\mu$ g/ml



For full product information, images and publications, please visit our [website](#).