

RIMS3 blocking peptide

Cat. No. GTX31334-PEP
Applications Neutralizing/Inhibition

Package
50 µg

PRODUCT

Summary RIMS3 peptide is used for blocking the activity of RIMS3 antibody (catalog number GTX31334).

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

Neutralizing/Inhibition	Assay dependent
-------------------------	-----------------

Not tested in other applications.

Properties

Form	Liquid
-------------	--------

Buffer	PBS, 0.1% BSA
---------------	---------------

Preservative	0.02% Sodium azide
---------------------	--------------------

Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
----------------	--

Concentration	200 µg/ml (Please refer to the vial label for the specific concentration.)
----------------------	--

Conjugation	Unconjugated
--------------------	--------------

Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>
-------------	---



For full product information, images and publications, please visit our [website](#).