

KDM3B antibody

Cat. No. GTX31386

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	2.5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	192 kDa. (Note)
Product Note	This antibody will not cross-react with JMJD1A or JMJD1C.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	KDM3B antibody was raised against a 20 amino acid synthetic peptide from near the center of human KDM3B. The immunogen is located within amino acids 1050 - 1100 of KDM3B.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.


DATA IMAGES



GTX31386 WB Image

WB analysis of rat liver tissue lysate using GTX31386 KDM3B antibody.


Dilution : 1 µg/ml (Lane A) and 2 µg/ml (Lane B)



GTX31386 IHC-P Image

IHC-P analysis of rat liver tissue using GTX31386 KDM3B antibody.

Working concentration : 20 µg/ml



GTX31386 IHC-P Image

IHC-P analysis of rat liver tissue using GTX31386 KDM3B antibody.

Working concentration : 2.5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).