

LMX1A blocking peptide

Cat. No. GTX31507-PEP

Applications	Neutralizing/Inhibition
--------------	-------------------------

References (1)

Package

50 µg

PRODUCT

Summary	LMX1A peptide is used for blocking the activity of LMX1A antibody (catalog number GTX31507).
---------	--

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
--------------------	----------------------

Neutralizing/Inhibition	Assay dependent
-------------------------	-----------------

Not tested in other applications.

Properties

Form	Liquid
------	--------

Buffer	PBS, 0.1% BSA
--------	---------------

Preservative	0.02% Sodium azide
--------------	--------------------

Storage	Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.
---------	--

Concentration	200 µg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Conjugation	Unconjugated
-------------	--------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
------	---



For full product information, images and publications, please visit our [website](#).