

GRIP1 antibody

Cat. No. GTX31521

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 122 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rabbit polyclonal GRIP1 antibody was raised against a 17 amino acid peptide near the carboxy terminus of human GRIP1. The immunogen is located within the last 50 amino acids of GRIP1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX31521 WB Image


WB analysis of 293 cell lysate in (A) the absence and (B) the presence of blocking peptide using GTX31521 GRIP1 antibody.

Working concentration : 1 µg/ml


GTX31521 IHC-P Image

IHC-P analysis of rat brain tissue using GTX31521 GRIP1 antibody.

Working concentration : 20 µg/ml


GTX31521 IHC-P Image

IHC-P analysis of rat brain tissue using GTX31521 GRIP1 antibody.

Working concentration : 2.5 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).