

RC3H1 antibody

Cat. No. GTX31528

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 126 kDa. ([Note](#))

Product Note At least three isoforms of RC3H1 are known to exist; this antibody will detect all three isoforms.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rabbit polyclonal RC3H1 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human RC3H1. The immunogen is located within amino acids 1040 - 1090 of RC3H1.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTx31528 WB Image

WB analysis of HeLa cell lysate using GTx31528 RC3H1 antibody.
Working concentration : 1 µg/ml



GTx31528 IHC-P Image

IHC-P analysis of human small intestine tissue using GTx31528 RC3H1 antibody.
Working concentration : 20 µg/ml



GTx31528 IHC-P Image

IHC-P analysis of human small intestine tissue using GTx31528 RC3H1 antibody.
Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](#).