

PRKCDBP antibody

Cat. No. GTX31530

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	28 kDa. (Note)
Product Note	It is predicted to not cross-react with other members of the cavin family.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Rabbit polyclonal PRKCDBP antibody was raised against an 17 amino acid peptide near the center of human PRKCDBP. The immunogen is located within amino acids 150 - 200 of PRKCDBP.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX31530 IHC-P Image


IHC-P analysis of human spleen tissue using GTX31530 PRKCDBP antibody.

Working concentration : 2.5 µg/ml


GTX31530 IHC-P Image

IHC-P analysis of human spleen tissue using GTX31530 PRKCDBP antibody.

Working concentration : 20 µg/ml


GTX31530 WB Image

WB analysis of A-20 cell lysate in (A) the absence and (B) the presence of blocking peptide using GTX31530 PRKCDBP antibody.

Working concentration : 1 µg/ml



For full product information, images and publications, please visit our [website](#).