

GBP5 antibody

Cat. No. GTX31537

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 µg/mL
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW	67 kDa. (Note)
Product Note	GBP5 antibody is predicted to not cross-react with other GBP family members.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	GBP5 antibody was raised against a 16 amino acid peptide near the carboxy terminus of human GBP5 .The immunogen is located within the last 50 amino acids of GBP5.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx31537 WB Image


WB analysis of rat liver tissue lysate in (A) the absence and (B) the presence of blocking peptide using GTx31537 GBP5 antibody.

Working concentration : 1 µg/ml


GTx31537 IHC-P Image

IHC-P analysis of human liver tissue using GTx31537 GBP5 antibody.

Working concentration : 2.5 µg/ml


GTx31537 IHC-P Image

IHC-P analysis of human liver tissue using GTx31537 GBP5 antibody.

Working concentration : 20 µg/ml



For full product information, images and publications, please visit our [website](#).