

GABARAP antibody

Cat. No. GTX31541

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human, Mouse, Rat

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	5 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 14 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	GABARAP antibody was raised against a 19 amino acid peptide near the amino terminus of human GABARAP. The immunogen is located within the first 50 amino acids of GABARAP.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX31541 IHC-P Image


IHC-P analysis of mouse brain tissue using GTX31541 GABARAP antibody.

Working concentration : 20 µg/ml

Green : Primary antibody


Blue : DAPI

Red : Actin


GTX31541 IHC-P Image


IHC-P analysis of mouse brain tissue using GTX31541 GABARAP antibody.

Working concentration : 5 µg/ml


GTX31541 WB Image

WB analysis of HeLa cell lysate using GTX31541 GABARAP antibody.

Working concentration : 1 µg/ml


GTX31541 IHC-P Image

IHC-P analysis of human brain tissue using GTX31541 GABARAP antibody.

Working concentration : 20 µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).