

IL1F6 antibody

Cat. No. GTX31555

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P, ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1 - 2 µg/mL
IHC-P	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 18 kDa. ([Note](#))

Product Note IL-36A antibody is human specific. IL-36A antibody will not cross-react with IL-36B or IL-36G.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	IL1F6 antibody was raised against an 11 amino acid peptide near the carboxy terminus of human IL1F6. The immunogen is located within the last 50 amino acids of IL1F6.
Purification	Purified by antigen-affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx31555 IHC-P Image

IHC-P analysis of human brain tissue using GTx31555 IL1F6 antibody.
Working concentration : 20 µg/ml


GTx31555 WB Image

WB analysis of MCF7 cell lysate in (A) the absence and (B) the presence of blocking peptide using GTx31555 IL1F6 antibody.
Working concentration : 1 µg/ml


GTx31555 IHC-P Image

IHC-P analysis of human brain tissue using GTx31555 IL1F6 antibody.
Working concentration : 5 µg/ml



For full product information, images and publications, please visit our [website](#).