

OLIG2 blocking peptide

Cat. No. GTX31569-PEP
Applications

Neutralizing /Inhibition

Package

50 µg

PRODUCT
Summary

OLIG2 peptide is used for blocking the activity of OLIG2 antibody (catalog number GTX31569). It usually blocks the antibody activity completely in Western blot by incubating the peptide with equal volume of antibody for 30 min at 37°C.

Applications
Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution
Recommended dilution

Neutralizing /Inhibition

Assay dependent

Not tested in other applications.

Properties

Form Liquid

Buffer PBS, 0.1% BSA

Preservative 0.02% Sodium azide

Storage Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration 200 µg/ml (Please refer to the vial label for the specific concentration.)

Conjugation Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).